 IMMULITE<sup>®</sup> 2000

**AFP**

For use on IMMULITE<sup>®</sup> 2000 systems

**SIEMENS**

## Symbol Definition

**CONTROL -**

En: Negative Control  
De: Negativkontrolle  
Es: Control Negativo  
Fr: Contrôle négatif  
It: Controllo negativo  
Pt: Controllo Negativo

**CONTROL AB**

En: Control Antibody  
De: Kontroll-Antikörper  
Es: Anticuerpo Control  
Fr: Anticorps du contrôle  
It: Anticorpo di controllo  
Pt: Anticorpo do Controllo

**PRE A**

En: Pretreatment  
Solution

**PRE B**

De: Vorbehandlungs-  
lösung  
Es: Solución de  
Pretratamiento  
Fr: Solution de  
prétraitement  
It: Soluzione di  
pretrattamento  
Pt: Solução de Pré-  
tratamento

**DITHIOTHREITOL**

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determinations of AFP are particularly valuable. The presence of residual tumor is strongly suggested if post-operative AFP concentrations fail to return to normal.<sup>9,15,23,29</sup> The accurate interpretation of post-surgical changes in AFP concentration requires consideration of its metabolic decay rate.<sup>21,22,24,25</sup> When utilizing AFP for monitoring therapy or disease recurrence during chemotherapy, it should be noted that levels often fall rapidly during chemotherapy, reaching normal levels while tumor masses are still evident.<sup>17,21</sup> In such instances, completion of the planned therapy has been recommended.<sup>21</sup>

Following therapy or surgery, serial measurements of AFP have also proved clinically useful when monitoring for progression or recurrence of disease in patients with nonseminomatous testicular cancer. It has been reported that AFP levels frequently rise during disease progression and fall during disease remission.<sup>9,17,19</sup> Elevated AFP levels have frequently been observed to accompany tumor recurrence before progressive disease is clinically evident.<sup>9,16</sup>

### Fetal Open Neural Tube Defects

AFP is detectable not only in fetal serum, but also in amniotic fluid and maternal serum. A concentration gradient exists such that when the fetal serum AFP level is 2000 kIU/mL, the amniotic fluid AFP (AFAP) level is 20 kIU/mL, and the maternal serum AFP (MSAFP) is 0.02 kIU/mL. In normal pregnancy, the fetal serum AFP concentration peaks at 14 weeks gestation.<sup>34</sup> The AFAP concentration peaks at about 12 weeks and the MSAFP peaks at approximately 28–32 weeks gestation.<sup>36</sup> The fall in AFAP concentration reflects the fall in fetal serum AFP concentration which results from increased fetal size and fluid volume.<sup>34</sup> Elevated levels of MSAFP and AFAP may occur most often due to multiple pregnancy and due to incorrect gestational age.

Measurement of AFP concentrations is clinically valuable in screening for open NTDs and other fetal abnormalities.<sup>35</sup> Pregnancies associated with open NTDs present with elevated AFP levels. Excess AFP gains access to amniotic fluid, and to a lesser extent to the maternal serum, by

transudation across the exposed surface of the fetus or across damaged glomeruli.<sup>35,37</sup> These conditions are found in open NTDs including open spina bifida and anencephaly, omphalocele, and congenital nephrosis.<sup>32,39</sup> Additional causes of elevated AFP concentrations including both maternal and fetal sources are impending spontaneous abortion, fetal distress or death, oligohydramnios, toxemia, gastroschisis, Meckel's syndrome, sacrococcygeal teratoma, Turner's syndrome and maternal hepatic and oncologic disorders.<sup>35</sup>

Recommended protocols for open NTD screening have been published.<sup>33,35</sup> The cutoff levels for maternal serum and amniotic fluid can be chosen to optimize the needs of the populations being tested based upon varying prevalence of open NTDs. Cutoffs commonly utilize multiples of the median of 2.0 or 2.5 for MSAFP and AFAP testing. The optimal time for screening MSAFP is between the 16th and 18th weeks of pregnancy, although screening is still effective before or after this period. Elevated AFP concentrations may be subjected to a repeat sampling and analysis to exclude transient rises.

More commonly, ultrasonography is employed to rule out multiple pregnancies and to confirm gestational age. Ultrasonography may also identify signs of open NTDs, particularly anencephaly which is a large, easily visualized lesion. If correction for gestational age or multiple pregnancy does not result in an AFP concentration within the normal range, then diagnostic ultrasonography and/or amniotic fluid sampling is indicated. The greatest diagnostic power can be achieved by combining biochemical analysis of amniotic fluid and diagnostic ultrasonography in cases of a positive MSAFP screen.<sup>35</sup>

Elevated MSAFP results are not diagnostic for NTDs and should not be considered a cause for termination of pregnancy. An overlap exists in the distributions of AFP concentrations from pregnancies with and without open NTDs. Closed NTDs, for example, are not usually associated with increased MSAFP or AFAP concentrations. Thus, further testing is required to define fetal status. In light of these considerations and the multiple causes for elevated AFP

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**Amniotic Fluid Dilution Factor: 100**  
 All amniotic fluid samples must first be diluted 1-in-101 using on-board Multi-Diluent 2 before being assayed. Select 100 in the Dilution Factor window.

**Storage**

**Serum:** 3 days at 2-8°C. Freeze at -20°C if not assayed within 3 days.

**Amniotic Fluid:** Amniotic fluid samples should be stored at -20°C. Aliquot if necessary to avoid repeated freezing and thawing. Allow the sample to come to room temperature (15-28°C) before assay, and mix by *gentle* swirling or inversion. Do not attempt to thaw specimens by heating them in a waterbath. If specimens are to be mailed, samples should be packed in dry ice if the time in transit exceeds 72 hours, or if elevated temperatures are a concern, as in warm climates or during the summer. If a repeat analysis is required, the original type of specimen should be taken to maintain consistency of results.

**Warnings and Precautions**

For *in vitro* diagnostic use.

**Reagents:** Store at 2-8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight (see insert).

**Water:** Use distilled or deionized water.

**Materials Supplied**

Components are a matched set. Labels on the inside box are needed for the assay.

**AFP Bead Pack (L2AP12)**

With barcode. 200 beads, coated with murine monoclonal anti-AFP. Stable at 2-8°C until expiration date.  
**L2KAP2:** 1 pack **L2KAP6:** 3 packs

**AFP Reagent Wedge (L2APAZ)**

With barcode. 11.5 mL of a protein buffer/nonhuman serum matrix; and 11.5 mL of alkaline phosphatase (bovine calf intestine) conjugated to polyclonal rabbit anti-AFP, in buffer. Stable at 2-8°C until expiration date.  
**L2KAP2:** 1 wedge **L2KAP6:** 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

**AFP Adjustors (L2APJ3, L2APJ4)**

Two vials (Low and High), 2.0 mL each, of AFP in a bovine serum matrix. Stable at 2-8°C for 30 days after opening, or for 6 months (aliquoted) at -20°C.  
**L2KAP2:** 1 set **L2KAP6:** 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes, so that the barcodes can be read by the on-board reader.

**Kit Components Supplied Separately****Multi-Diluent 2 (L2M2Z, L2M2Z4)**

For the on-board dilution of high serum samples and for amniotic fluid samples. One vial of concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Stable at 2-8°C for 30 days after opening, or for 6 months (aliquoted) at -20°C.  
**L2M2Z:** 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 x 100 mm test tube, so that the barcode can be read by the on-board reader.  
**L2M2Z:** 3 labels **L2M2Z4:** 5 labels

**Analysis of amniotic fluid requires a 1-in-101 dilution of the sample (on-board dilution with Multi-Diluent 2).**

**L2SUBM:** Chemiluminescent Substrate  
**L2PWSM:** Probe Wash  
**L2KPM:** Probe Cleaning Kit

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serum AFP levels. However, elevated AFP levels have been observed in patients diagnosed with seminomas accompanied by metastases of nonseminomatous testicular cancer.<sup>8</sup>

A significant increase of AFP levels in patients considered free of metastatic tumor may indicate the development of metastasis. Elevated levels after surgery may indicate incomplete removal of the tumor or the presence of metastases.

Elevated levels of serum AFP are associated with benign liver conditions such as hepatitis and cirrhosis. Most (95%) of the patients with these benign diseases have AFP levels lower than 200 ng/mL (165 IU/mL).<sup>8-15</sup>

### AFP Values in Maternal Serum and Amniotic Fluid

Due to potential variation in testing at different laboratories, it is recommended that a particular testing center determine its own set of median AFP values for weeks 15 to 20 of gestation, measured in the population to be screened. Cutoff values commonly utilize multiples of the medians (MoM) of 2.0 or 2.5 for maternal serum and amniotic fluid testing. Each AFP test result can then be expressed as a multiple of the unaffected population median value. This is obtained by dividing the AFP value by the median value for its corresponding gestational week. Gestational weeks are defined as completed gestational weeks; e.g., 16 weeks, 6 days would be considered the 16th week. It has been recommended that median and MoM values determined for each gestational week be based upon at least 100 maternal sera and 50 amniotic fluids from unaffected singleton pregnancies with confirmed gestational age.

Provided below are medians for maternal serum samples, calculated by a weighted log-linear regression from data collected from unaffected, singleton pregnancies at three clinical sites in the United States:

Gestational Week	No. of Specimens	Medians IU/mL*	Multiples of Regressed Medians (IU/mL)		
			2.0	2.5	3.0
15	370	24.9	49.8	62.3	74.7
16	605	28.5	57.0	71.3	85.5
17	589	32.6	65.2	81.5	97.8
18	431	37.2	74.4	93.0	111.6
19	221	42.5	85.0	106.3	127.5
20	91	48.6	97.2	121.5	145.8

\*Regressed

Provided below are medians for amniotic fluid samples, calculated by a weighted log-linear regression from data collected from unaffected, singleton pregnancies at two clinical sites in the United States:

Gestational Week	No. of Specimens	Medians kIU/mL*	Multiples of Regressed Medians (kIU/mL)		
			2.0	2.5	3.0
15	76	13.0	26.0	32.5	39.0
16	89	10.7	21.4	26.8	32.1
17	53	8.73	17.5	21.8	26.2
18	54	7.14	14.3	17.9	21.4
19	46	5.84	11.7	14.6	17.5
20	23	4.78	9.56	12.0	14.3

\*Regressed

### Limitations

**Diagnosis:** The occurrence of elevated serum AFP levels in conditions other than nonseminomatous testicular cancer precludes the use of AFP measurements in the diagnosis of nonseminomatous testicular cancer.

**Screening:** AFP measurements can not be recommended as a screening procedure to detect cancer in the general population. Elevated concentrations of serum AFP have been observed not only in patients with nonseminomatous testicular cancer but also in malignant conditions such as hepatocellular carcinoma, ovarian cancer, and gastrointestinal and pulmonary cancer. Benign hepatic conditions such as acute viral hepatitis, chronic active hepatitis and cirrhosis may present with elevated concentrations of serum AFP. Elevated AFP concentrations have also been observed in pregnancy, ataxia telangiectasia and hereditary tyrosinemia.

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IMMULITE 2000 AFP (PIL2KAP-21, 2014-02-26)

The assay was also compared to IMMULITE AFP on maternal serum samples, in the range from approximately 10 to 120 IU/mL. By linear regression: (IML 2000) = 1.01 (IML) + 0.154 IU/mL  
r = 0.982  
n = 346

Means:  
33.8 IU/mL (IML)  
34.3 IU/mL (IML 2000)

95% Confidence Interval (CI)	Slope	Intercept
Lower	0.89	0.94
Upper	0.93	2.68

(IML 2000) = 0.73 (Kit B) + 5.22 IU/mL  
r = 0.97  
n = 1015

95% Confidence Interval (CI)	Slope	Intercept
Lower	0.72	4.66
Upper	0.74	5.79

In one of the the studies above, IMMULITE 2000 AFP results were compared to Kit B in a linear regression for amniotic fluid samples, in the range from nondetectable to 286 kIU/mL\*. (IML 2000) = 0.79 (Kit B) + 2.27 kIU/mL  
r = 0.99  
n = 200

95% Confidence Interval (CI)	Slope	Intercept
Lower	0.77	1.78
Upper	0.81	2.76

\* Amniotic fluid samples were diluted 1-in-101 automatically by the IMMULITE 2000 instrument.

The assay was also compared to IMMULITE AFP on amniotic fluid samples, in the range from approximately 3 to 20 kIU/mL\*. (See graph.) By linear regression:

(IML 2000) = 1.03 (IML) + 0.52 kIU/mL  
r = 0.96  
n = 46

Means:  
10.0 kIU/mL (IML)  
10.8 kIU/mL (IML 2000)

95% Confidence Interval (CI)	Slope	Intercept
Lower	0.93	-0.50
Upper	1.12	1.54

\* Amniotic fluid samples were diluted 1-in-101 automatically by the IMMULITE 2000 instrument.

### Clinical Sensitivity for Maternal Serum, n = 9:

Gestational Week	% > 2.0 MoM	% > 2.5 MoM	% > 3.0 MoM
15-20	100%	77.8%	66.7%
95% CI for All Samples	66.4%–100%	40.0%–97.2%	29.9%–92.5%

### Clinical Specificity for Maternal Serum:

Gestational Week	n	% ≤ 2.0 MoM	% ≤ 2.5 MoM	% ≤ 3.0 MoM
15	276	94.2%	97.5%	98.6%
16	304	96.1%	99.0%	99.7%
17	272	97.1%	99.3%	99.6%
18	287	95.8%	98.6%	99.3%
19	152	93.4%	98.0%	99.3%
20	41	95.1%	100%	100%
15-20	1332	95.5%	98.6%	99.3%
95% CI for All Samples		94.2%–96.5%	97.8%–99.1%	98.7%–99.7%

### Clinical Sensitivity for Amniotic Fluid, n = 8:

Gestational Week	% > 2.0 MoM	% > 2.5 MoM	% > 3.0 MoM
15-20	87.5%	87.5%	87.5%
95% CI for All Samples	47.3%–99.7%	47.3%–99.7%	47.3%–99.7%

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