Free Beta HCG

For use on IMMULITE® 2000 systems

SIEMENS

Symbol Definition
DIL
 En: Sample Diluent
 De: Probenverdünnungsreagenz
 Es: Diluyente para muestras
 Fr: Diluant échantillon
 It: Diluente per Campioni
 Pt: Diluente de Amostra

Symbol Definition
CONTROL +
 En: Control
 De: Kontrolle
 Es: Control
 Fr: Contrôle
 It: Controllo
 Pt: Controllo

Symbol Definition
CONTROL 1
CONTROL 2
CONTROL 3

Symbol Definition
CONTROL -
 En: Negative Control
 De: Negativkontrolle
 Es: Control Negativo
 Fr: Contrôle négatif
 It: Controllo negativo
 Pt: Controllo Negativo

Symbol Definition
CONTROL AB
 En: Control Antibody
 De: Kontroll-Antikörper
 Es: Anticuerpo Control
 Fr: Anticorps du contrôle
 It: Anticorpo di Controllo
 Pt: Anticorpo do Controllo

Symbol Definition
PRE A
PRE B

Symbol Definition
DITHIOTHREITOL
 En: Dithiothreitol
 Solution
 De: Dithiothreitol-Lösung
 Es: Solución de Ditiotretol
 Fr: Solution de Dithiothreitol
 It: Soluzione di Ditiotretolo
 Pt: Solução de Ditiotretol

Symbol Definition
BORATE-KCN BUF
 En: Borate-KCN Buffer Solution
 De: Borat-KCN-Puffer
 Es: Solución Tampón Borato-KCN
 Fr: Solution tampon Borate-Cyanure de Potassium
 It: Soluzione Tampone Borato-KCN
 Pt: Solução Tamponzada de Borato-KCN

chemiluminescent substrate is added to the reaction tube and the signal is generated in proportion to the bound enzyme.

Incubation Cycles: 2 x 30 minutes.
Time to First Result: 65 minutes.

Specimen Collection
 The use of an ultracentrifuge is recommended to clear lipemic samples. Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Free Beta HCG has not been tested with all possible variations of tube types.

Volume Required: 5 µL serum.
Storage: 7 days at 2–8°C or longer at –20°C. Avoid repeated freeze/thaw cycles.

Dilution of High Samples: Samples expected to have free β-HCG levels above the assay's reportable range should be diluted with HCG Sample Diluent (L2CGZ) before assay.

Warnings and Precautions
 For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD
 Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2), as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412
P273, P501
 Harmful to aquatic life with long lasting effects.
 Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations.
Contains: sodium azide; Free Beta HCG Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.
 Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)
Water: Use distilled or deionized water.

Specificity: The antibody is highly specific for the free β subunit of HCG. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/dL has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 600 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 22 volunteers into serum, lithium-heparin, EDTA and Becton Dickinson SST[®] plastic vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of free β -HCG to obtain values throughout the reportable range of the assay, and then assayed by the IMMULITE 2000 Free β -HCG procedure.

$r = 0.985, n = 22$
(LI Heparin) = 0.88 (Serum) + 3.05 ng/mL
(EDTA) = 0.89 (Serum) + 2.02 ng/mL
 $r = 0.983, n = 22$
(SST) = 0.99 (Serum) + 0.95 ng/mL
 $r = 0.982, n = 22$

Means:
91.1 ng/mL (Serum)
83.3 ng/mL (LI Heparin)
82.8 ng/mL (EDTA)
91.1 ng/mL (SST)

Method Comparison: The assay was compared to a commercially available Free Beta HCG assay (Kryptor) on 395 samples. (Concentration range: approximately 5.7 to 197 ng/mL. See graph.) By linear regression:
(IML 2000) = 1.10 (Kryptor) - 0.25 ng/mL
 $r = 0.985$

Means:
45.8 ng/mL (IMMULITE 2000)
41.7 ng/mL (Kryptor)

and testicular tumors. Thus, the assay is not applicable for this clinical utility.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference, however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data
See Tables and Graphs for data representative of the assay's performance. Results are expressed in ng/mL, where 1 ng/mL of free β -HCG is equivalent to 1 mIU/mL of free β -HCG in terms of the World Health Organization's First International Reference Preparation of Chorionic Gonadotrophin Beta Subunit Human, number 75/551 (1st IRP 75/551). Unless otherwise noted, all results were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.

Reportable Range: 2 to 200 ng/mL
Analytical Sensitivity: 1 ng/mL

High-dose Hook Effect: None up to at least 133,735 ng/mL

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table.)

Recovery: Samples spiked 1 to 19 with three free β -HCG solutions (50, 500 and 1,500 ng/mL) were assayed. (See "Recovery" table.)

Symbol Definition
En: Upper limit of temperature ($\leq -20^{\circ}\text{C}$)
De: Obere Temperaturgrenze ($\leq -20^{\circ}\text{C}$)
Es: Limite superior de temperatura ($\leq -20^{\circ}\text{C}$)
Fr: Limite supérieure de température ($\leq -20^{\circ}\text{C}$)
It: Limite superiore di temperatura ($\leq -20^{\circ}\text{C}$)
Pt: Limite máximo de temperatura ($\leq -20^{\circ}\text{C}$)

En: Lower limit of temperature ($\geq 2^{\circ}\text{C}$)
De: Mindesttemperatur ($\geq 2^{\circ}\text{C}$)
Es: Limite inferior de temperatura ($\geq 2^{\circ}\text{C}$)
Fr: Limite inférieure de température ($\geq 2^{\circ}\text{C}$)
It: Limite inferiore di temperatura ($\geq 2^{\circ}\text{C}$)
Pt: Limite mínimo de temperatura ($\geq 2^{\circ}\text{C}$)

En: Do not freeze ($> 0^{\circ}\text{C}$)
De: Nicht einfrieren ($> 0^{\circ}\text{C}$)
Es: No congelar ($> 0^{\circ}\text{C}$)
Fr: Ne pas congeler ($> 0^{\circ}\text{C}$)
It: Non congelare ($> 0^{\circ}\text{C}$)
Pt: Não congelar ($> 0^{\circ}\text{C}$)

En: Do not reuse
De: Nicht zur Wiederverwendung
Es: No reutilizar
Fr: Ne pas réutiliser
It: Non riutilizzare
Pt: Não reutilizar

En: Keep away from sunlight
De: Vor Sonneneinstrahlung schützen
Es: Proteger de la luz solar
Fr: Maintenir hors de portée de la lumière du soleil
It: Non esporre alla luce del sole
Pt: Manter afastado da luz solar

En: Batch code
De: Chargenbezeichnung
Es: Número de lote
Fr: Numéro de code du lot
It: Codice lotto
Pt: Código de lote

Symbol Definition
De: CE Mark
Es: Marca CE
Fr: Marque CE
It: Marchio CE
Pt: Marca CE

En: CE Mark with identification number of notified body
De: CE-Kennzeichen mit Identifikationsnummer der benannten Stelle
Es: Marca CE con número de identificación del organismo notificado
Fr: Marque CE avec numéro d'identification du corps notifié
It: Marchio CE con numero identificativo dell'ente notifico
Pt: Marca CE, com número de identificação do organismo notificado

En: Consult instructions for use
De: Bedienungsanweisung beachten
Es: Consulte las instrucciones de uso
Fr: Consulter le mode d'emploi
It: Consultare le istruzioni per l'uso
Pt: Consulte as instruções de utilização

En: Caution! Potential Biohazard
De: Vorsicht! Biologisches Risikomaterial
Es: ¡Precaución! Riesgo biológico potencial
Fr: Attention! Potentiel biologique dangereux
It: Attenzione! Potenziale Pericolo Biologico
Pt: Atenção! Potenciais Riscos Biológicos

En: Temperature limitation (2-8°C)
De: Temperaturgrenze (2-8°C)
Es: Limitación de temperatura (2-8°C)
Fr: Limites de température (2-8°C)
It: Limiti di temperatura (2-8°C)
Pt: Limites de temperatura (2-8°C)

Symbol Definition
En: Upper limit of temperature ($\leq -20^{\circ}\text{C}$)
De: Obere Temperaturgrenze ($\leq -20^{\circ}\text{C}$)
Es: Limite superior de temperatura ($\leq -20^{\circ}\text{C}$)
Fr: Limite supérieure de température ($\leq -20^{\circ}\text{C}$)
It: Limite superiore di temperatura ($\leq -20^{\circ}\text{C}$)
Pt: Limite máximo de temperatura ($\leq -20^{\circ}\text{C}$)

En: Lower limit of temperature ($\geq 2^{\circ}\text{C}$)
De: Mindesttemperatur ($\geq 2^{\circ}\text{C}$)
Es: Limite inferior de temperatura ($\geq 2^{\circ}\text{C}$)
Fr: Limite inférieure de température ($\geq 2^{\circ}\text{C}$)
It: Limite inferiore di temperatura ($\geq 2^{\circ}\text{C}$)
Pt: Limite mínimo de temperatura ($\geq 2^{\circ}\text{C}$)

En: Do not freeze ($> 0^{\circ}\text{C}$)
De: Nicht einfrieren ($> 0^{\circ}\text{C}$)
Es: No congelar ($> 0^{\circ}\text{C}$)
Fr: Ne pas congeler ($> 0^{\circ}\text{C}$)
It: Non congelare ($> 0^{\circ}\text{C}$)
Pt: Não congelar ($> 0^{\circ}\text{C}$)

En: Do not reuse
De: Nicht zur Wiederverwendung
Es: No reutilizar
Fr: Ne pas réutiliser
It: Non riutilizzare
Pt: Não reutilizar

En: Keep away from sunlight
De: Vor Sonneneinstrahlung schützen
Es: Proteger de la luz solar
Fr: Maintenir hors de portée de la lumière du soleil
It: Non esporre alla luce del sole
Pt: Manter afastado da luz solar

En: Batch code
De: Chargenbezeichnung
Es: Número de lote
Fr: Numéro de code du lot
It: Codice lotto
Pt: Código de lote

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El Sistema de Calidad de Siemens Healthcare Diagnostics Products Ltd. está certificado por la ISO 13485:2003.