 IMMULITE®
2000

HCG

For use on IMMULITE® 2000 systems

SIEMENS

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.^{23,25}

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412
P273, P501

Harmful to aquatic life with long lasting effects.
Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations.
Contains: sodium azide; HCG Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

HCG Bead Pack (L2CG12)

With barcode, 200 beads, coated with monoclonal murine anti-HCG.
Stable at 2–8°C until expiration date.
L2KCG2: 1 pack
L2KCG6: 3 packs

HCG Reagent Wedge (L2CGA2)

With barcode, 11.5 mL of alkaline phosphatase (bovine calf intestine) conjugated to polyclonal ovine anti-HCG, in buffer, with preservative. Stable at 2–8°C until expiration date.
L2KCG2: 1 wedge
L2KCG6: 3 wedges

Before use, tear off the top of the label at the perforation, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

HCG Adjustors (L2GL, L2GH)

Two vials (Low and High) 2 mL each, of HCG in an HCG-free human serum matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.
L2KCG2: 1 set
L2KCG6: 2 sets

Before running adjustors, place the appropriate Aliquot Labels (supplied with the kit) on test tubes, so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

HCG Sample Diluent (L2CGZ, L2CGZ4)

For the on-board dilution of high samples. One vial of concentrated (ready-to-use) HCG-free human serum, with barcode labels. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

Before use, place an appropriate label on a 16 x 100 mm test tube, so that the barcodes can be read by the on-board reader.
L2CGZ, L2CGZ4: 5 labels
L2CGZ: 50 mL
L2CGZ4: 100 mL

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash
L2KPM: Probe Cleaning Kit

Symbol Definition
 En: Use by
 De: Verwendbar bis
 Es: Fecha de caducidad
 Fr: A utiliser avant
 It: Usare entro
 Pt: Usar até

En: Health Hazard
 De: Gesundheitsgefährdung
 Es: Peligro para la salud
 Fr: Dangereux pour la santé
 It: Pericolo per la salute
 Pt: Perigo para a saúde

En: Exclamation Mark
 De: Ausrufezeichen
 Es: Signo de exclamación
 Fr: Point d'exclamation
 It: Punto esclamativo
 Pt: Ponto de exclamação

En: Corrosion
 De: Korrosion
 Es: Corrosion
 Fr: Corrosion
 It: Corrosione
 Pt: Corrosão

En: Skull and Crossbones
 De: Totenkopf mit gekreuzten Knochen
 Es: Calavera y tibias cruzadas
 Fr: Tête de mort sur tibias croisées
 It: Teschio e tibie incrociate
 Pt: Caveira sobre tibias cruzadas

En: Environment
 De: Umwelt
 Es: Medio ambiente
 Fr: Environnement
 It: Ambiente
 Pt: Ambiente

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CONTROL 3

En: Control

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Pathological Conditions

Ectopic pregnancies or abnormal intrauterine pregnancies have circulating HCG levels below normal for that gestational period; after an initial rise these levels plateau or fall.^{3,17,18}

Limitations

The procedure for urine samples is designed solely as a screen for pregnancy, using a 30 mIU/mL cutoff. Hence urinary results should be reported in a strictly qualitative manner, e.g. as "positive" or "negative" for pregnancy — concentration values should *not* be reported for this sample type.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See tables and graphs for data *representative* of the assay's performance. Results are expressed as mIU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Reportable Range: Up to 5000 mIU/mL, in terms of the WHO 3rd IS 75/537 for quantitative determinations in serum.

Analytical Sensitivity: 0.4 mIU/mL

High-Dose Hook Effect: None up to 600,000 mIU/mL

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three solutions (1313, 5720 and 30,340 mIU/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The assay is highly specific for HCG, with low cross-reactivity to other glycoprotein hormones present in patient samples. (See "Specificity" table.)

For concentrations as high as 10,000 ng/mL, the HCG α -subunit (WHO 75/569) showed less than 1% cross-reactivity. HCG β -subunit (WHO 75/551), spiked into HCG Sample Diluent, yielded the tabulated results. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/dL has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 384 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of lipemia in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 14 volunteers into plain, heparinized, EDTA vacutainer tubes, and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of HCG, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 HCG procedure.

(Heparin) = 1.06 (Serum) – 7.06 mIU/mL
 $r = 0.99$
 (EDTA) = 1.08 (Serum) – 6.7 mIU/mL
 $r = 0.99$
 (SST) = 0.98 (PlainTubes) – 3.25 mIU/mL
 $r = 0.99$

IMMULITE is a trademark of Siemens Healthcare Diagnostics. © 2015 Siemens Healthcare Diagnostics. All rights reserved.

Made in: UK

 Siemens Healthcare Diagnostics Products Ltd
Glyn Rhonwy, Llanberis,
Gwynedd LL55 4EL
United Kingdom

0088

When used for Trisomy 21 testing bei Verwendung für Tests auf Trisomie 21 cuando se usa para el análisis de Trisomía 21 en cas d'utilisation pour les tests de la trisomie 21 se utilizzato per i test della trisomia 21 quando utilizzato para testes de Trisomia 21

2015-05-21
PLI2KCG – 33

Changes in this Edition:

cf#EU2534: 1) Removed references to discontinued products. 2) Under Warnings and Precautions, added biohazard symbol and potential biohazard statement, caution statement, required GHS hazard and precautionary codes and statements, offering chemicals and sources. 3) Under References, added citations from the potential biohazard statement. 4) Changed "Origin: UK" to "Made in: UK". 5) Updated the manufacturer address to include Glyn Rhonwy. 6) In Understanding the Symbols, revised the symbol for catalog number, added a GHS pictogram for health hazard, and updated existing symbols with GHS pictograms and terminology for exclamation mark, corrosion, skull and crossbones, and environment.

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immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8. 24) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections: Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute. 2005. NCCLS Document M29-A3. 25) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard. 29 CFR 1910.1030.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor. www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

Tables and Graphs

Precision (mIU/mL)

	Within-Run ¹			Total ²	
	Mean ³	SD ⁴	CV ⁶	SD	CV
1	6.5	0.43	6.6%	0.48	7.4%
2	103	3.3	3.2%	4.6	4.5%
3	239	7.6	3.2%	11.1	4.6%
4	438	13	3.0%	21	4.8%
5	3120	79	2.5%	150	4.8%

Linearity (mIU/mL)

Dilution ¹	Observed ²	Expected ³	%OE ⁴
1	16 in 16 ⁵	204	—
	8 in 16	100	102
	4 in 16	50	51
	2 in 16	24	26
	1 in 16	12	13
2	16 in 16	552	—
	8 in 16	273	276
	4 in 16	134	138
	2 in 16	68	69
	1 in 16	33	35
3	16 in 16	1144	—
	8 in 16	556	572
	4 in 16	267	286
	2 in 16	140	143
	1 in 16	70	72
4	16 in 16	2088	—
	8 in 16	1020	1044
	4 in 16	503	522
	2 in 16	240	261
	1 in 16	119	131

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El Sistema de Calidad de Siemens Healthcare Diagnostics Products Ltd. está certificado por la ISO 13485:2003.