PAPP-A

For use on IMMULITE® 2000 systems

SIEMENS

Symbol Definition

CONTROL +

En: Positive Control
De: Positivkontrolle
Es: Control Positivo
Fr: Contrôle positif
It: Controllo positivo
Pt: Controllo Positivo

PRE A
PRE B

En: Low Positive Control
De: Schwachpositivkontrolle
Es: Control Positivo bajo
Fr: Contrôle positif faible
It: Controllo Positivo Basso
Pt: Controllo Positivo Baixo

CONTROL + L

En: Negative Control
De: Negativkontrolle
Es: Control Negativo
Fr: Contrôle négatif
It: Controllo negativo
Pt: Controllo Negativo

CONTROL -

En: Negative Control
De: Negativkontrolle
Es: Control Negativo
Fr: Contrôle négatif
It: Controllo negativo
Pt: Controllo Negativo

CONTROL AB

En: Control Antibody
De: Kontroll-Antikörper
Es: Anticuerpo Control
Fr: Anticorps du contrôle
It: Anticorpo di Controllo
Pt: Anticorpo do Controllo

Symbol Definition

PRE A
PRE B

En: Pretreatment Solution
De: Vorbehandlungslosung
Es: Solucion de Pretratamiento
Fr: Solution de prettraitement
It: Soluzione di prettrattamento
Pt: Solução de Pré-tratamento

DITHIOHREITOL

En: Dithiothreitol Solution
De: Dithiothreitol-Lösung
Es: Solucion de Ditiotretol
Fr: Solution de Dithiothreitol
It: Soluzione di Ditiotretolo
Pt: Solução de Ditiotretol

BORATE-KCN BUF

En: Borate-KCN Buffer Solution
De: Borat-KCN-Puffer
Es: Solucion Tampón Borato-KCN
Fr: Solution tampon Borate-Cyanure de Potassium
It: Soluzione Tampone Borato-KCN
Pt: Solução Tamponzada de Borato-KCN

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.⁸⁻¹⁰

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412
P273, P501

Harmful to aquatic life with long lasting effects.
Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations.
Contains: 2-methyl-2H-isothiazol-3-one, PAPP-A Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen, and for antibodies to hepatitis C.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

PAPP-A Bead Pack (L2PC12)

With barcode. 200 beads, coated with monoclonal murine anti-PAPP-A. Stable at 2–8°C until expiration date.
L2KPC2: 1 pack

PAPP-A Reagent Wedge (L2PCA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-PAPP-A in buffer. Stable at 2–8°C until expiration date.
L2KPC2: 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge, snap the sliding cover down into the ramps on the reagent lid.

PAPP-A Adjustors (LPCL, LPCH)

Two vials (Low and High) of lyophilized PAPP-A, in a nonhuman serum matrix. Reconstitute each vial with 2.0 mL distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 3 months (aliquotted) at –20°C.
L2KPC2: 1 set

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Multi-Diluent 2 (L2M2Z, L2M2Z4)

For the on-board dilution of patient samples. One vial, concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.
L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 x 100 mm test tube, so that the barcodes can be read by the on-board reader.
L2M2Z: 3 labels **L2M2Z4:** 5 labels

L2SUBM: Chemiluminescent Substrate
L2PWSM: Probe Wash
L2KPM: Probe Cleaning Kit
L2RXT: Reaction Tubes (disposable)

0088

Symbol Definition

En: CE Mark with identification number of notified body
De: CE-Kennzeichen mit Identifikationsnummer der benannten Stelle
Es: Marca CE con número de identificación del organismo notificado
Fr: Marque CE avec numéro d'identification du corps notifié
It: Marchio CE con numero identificativo dell'ente notificato
Pt: Marca CE, com número de identificação do organismo notificado

En: Consult instructions for use
De: Bedienungsanweisung beachten
Es: Consulte las instrucciones de uso
Fr: Consulter le mode d'emploi
It: Consultare le istruzioni per l'uso
Pt: Consulte as instruções de utilização

En: Caution! Potential Biohazard
De: Vorsicht! Biologisches Risikomaterial
Es: ¡Precaución! Riesgo biológico potencial
Fr: Avertissement ! Risque biologique potentiel
It: Attenzione! Potenziale Pericolo Biologico
Pt: Atenção! Potenciais Riscos Biologicos

En: Temperature limitation (2-8°C)
De: Temperaturgrenze (2-8°C)
Es: Limitación de temperatura (2-8°C)
Fr: Limites de température (2-8°C)
It: Limiti di temperatura (2-8°C)
Pt: Limites de temperatura (2-8°C)

En: Upper limit of temperature ($\leq -20^{\circ}\text{C}$)
De: Obere Temperaturgrenze ($\leq -20^{\circ}\text{C}$)
Es: Límite superior de temperatura ($\leq -20^{\circ}\text{C}$)
Fr: Limite supérieure de température ($\leq -20^{\circ}\text{C}$)
It: Limite superiore di temperatura ($\leq -20^{\circ}\text{C}$)
Pt: Limite máximo de temperatura ($\leq -20^{\circ}\text{C}$)

**Symbol Definition**

En: Lower limit of temperature ($\geq 22^{\circ}\text{C}$)
De: Mindesttemperatur ($\geq 22^{\circ}\text{C}$)
Es: Límite inferior de temperatura ($\geq 22^{\circ}\text{C}$)
Fr: Limite inférieure de température ($\geq 22^{\circ}\text{C}$)
It: Limite inferiore di temperatura ($\geq 22^{\circ}\text{C}$)
Pt: Limite mínimo de temperatura ($\geq 22^{\circ}\text{C}$)

En: Do not freeze ($> 0^{\circ}\text{C}$)
De: Nicht einfrieren ($> 0^{\circ}\text{C}$)
Es: No congelar ($> 0^{\circ}\text{C}$)
Fr: Ne pas congeler ($> 0^{\circ}\text{C}$)
It: Non congelare ($> 0^{\circ}\text{C}$)
Pt: Não congelar ($> 0^{\circ}\text{C}$)

En: Do not reuse
De: Nicht zur Wiederverwendung
Es: No reutilizar
Fr: Ne pas réutiliser
It: Non riutilizzare
Pt: Não reutilizar

En: Keep away from sunlight
De: Vor Sonneneinstrahlung schützen
Es: Proteger de la luz solar
Fr: Maintenir hors de portée de la lumière du soleil
It: Non esporre alla luce del sole
Pt: Manter afastado da luz solar

En: Batch code
De: Chargenbezeichnung
Es: Número de lote
Fr: Numéro de code du lot
It: Codice lotto
Pt: Código de lote

En: Contains sufficient for (n) tests
De: Es reicht für (n) Tests
Es: Contiene suficiente para (n) pruebas
Fr: Contient du matériel suffisant pour (n) tests
It: Contiene materiale sufficiente per (n) test
Pt: Contém o suficiente para (n) testes

2008-01

En: Date format (year-month)
De: Datumformat (Jahr-Monat)
Es: Formato de fecha (año-mes)
Fr: Format de la date (année-mois)
It: Formato data (anno-mese)
Pt: Formato de data (ano-mês)

range of the assay, and then assayed by the IMMULITE 2000 PAPP-A procedure. (Heparin) = 1.08 (Serum) + 0.08 mIU/mL $r = 0.996$

Means:
4.72 mIU/mL (Serum)
5.18 mIU/mL (Heparin)

In another study, blood was collected from 24 volunteers into plain and EDTA tubes. Equal volumes of the matched samples were spiked with various concentrations of PAPP-A, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 PAPP-A procedure.
(EDTA) = 1.21 (Serum) + 0.02 mIU/mL $r = 0.995$

Means:
3.65 mIU/mL (Serum)
4.43 mIU/mL (EDTA)

EDTA plasma is not suitable for use.

Method Comparison 1: The assay was compared to a commercially available assay for PAPP-A (Kit A) on 152 patient samples. (Concentration range: approximately 0.3 to 10 mIU/mL. See graph.) By linear regression:
(IML 2000) = 0.99 (Kit A) + 0.5 mIU/mL $r = 0.973$

Means:
3.82 mIU/mL (IMMULITE 2000)
3.36 mIU/mL (Kit A)

Method Comparison 2: The assay was compared to IMMULITE PAPP-A on 283 samples. (Concentration range: approximately 0.3 to 10 mIU/mL. See graph.) By linear regression:
(IML 2000) = 1.03 (IML) - 0.08 mIU/mL $r = 0.995$

Means:
3.70 mIU/mL (IMMULITE 2000)
3.64 mIU/mL (IMMULITE)

References

- Lin, T-M, *et al.* Characterization of four human pregnancy-associated plasma proteins. *Am J Obstet Gynecol* 1974; 118: 223-236; 2). Lin, T-M, *et al.* Measurement of pregnancy-associated plasma proteins during human gestation. *J Clin Invest* 1974; 54: 576-82; 3) Pederson, J. F. *et al.* Serum levels of human placental lactogen, pregnancy-associated plasma protein A and endometrial secretory PP14 in first trimester of diabetic pregnancy. *Acta Obstet Gynecol Scand* 1998; 77: 155-8; 4) Wald NJ, *et al.* Serum screening for Down's syndrome between 8 and

14 weeks of pregnancy. *Br J Obstet Gynaecol* 1996; 103: 407-12; 5) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4; Wayne, PA: NCCLS, 1996; 6) Bindra R, *et al.* One-stop clinic for assessment of risk for trisomy 21 at 11-14 weeks: a prospective study of 15 030 pregnancies. *Ultrasound Obstet Gynecol* 2002; 3: 219-25; 7) Wald NJ, Hackshaw AK. Advances in antenatal screening for Down syndrome. *Baillieres Best Pract Res Clin Obstet Gynecol* 2000; 14(4): 563-80; 8) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR* 1988; 37: 377-82; 367-8; 9) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.

10) Federal Occupational Safety and Health Administration. *Bloodborne Pathogens Standard*, 29 CFR 1910.1030.

Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

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IMMULITE is a trademark of
All rights reserved.
El Sistema de Calidad de Siemens Healthcare Diagnostics Products Ltd. está certificado por la ISO 13485:2003.