Unconjugated Estriol

For use on IMMULITE® 2000 systems

SIEMENS

Symbol Definition

En: Pretreatment
Solution
De: Vorbehandlungs-
lösung
Es: Solución de
Pretratamiento
Fr: Solution de
prétraitement
It: Soluzione di
pretrattamento
Pt: Solução de
Pré-tratamento

PRE A

PRE B

En: Dithiothreitol
Solution
De: Dithiothreitol-
Lösung
Es: Solución de
Ditiotreitól
Fr: Solution de
Dithiothreitol
It: Soluzione di
Ditiotreitolo
Pt: Solução de
Ditiotreitól

DITHIOTHREITOL

En: Borate-KCN
Buffer Solution
De: Borat-KCN-Puffer
Es: Solución Tampón
Borato-KCN
Fr: Solution tampon
Borate-Cyanure de
Potassium
It: Soluzione
Tampone Borato-KCN
Pt: Solução
Tamponzada de
Borato-KCN

BORATE-KCN BUF

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples. Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution. Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time. High patient samples can be manually diluted with an appropriate diluent. Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Unconjugated Estriol has not been tested with all possible variations of tube types.

Volume Required: 40 µL serum

Storage: 7 days at 2–8°C or 6 months at –20°C. Avoid multiple freeze thaw cycles.

Although serum samples have been reported to be stable at –20°C, unpublished observations suggest a small increase of approximately 10–15% may occur during storage at this temperature. Therefore it is preferable to transport and store samples at 2–8°C and assay within one week of collection.

It is important that all samples drawn for serial monitoring of a patient's unconjugated estriol level be measured using the same assay system, and that all be subjected to the same storage conditions.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.^{28,30}

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H318, H315, H400, H411

P280, P273, P305 + P351 + P310, P501

Danger! Causes serious eye damage. Causes skin irritation. Very toxic to aquatic life. Toxic to lasting effects.

Wear protective gloves/protective clothing/eye protection/face protection. Avoid release to the environment. IF IN EYES: Rinse cautiously with water for several minutes. Immediately call a

POISON CENTER or doctor/physician. Dispose of contents and container in accordance with all local, regional, and national regulations. Contains: zinc chloride; Unconjugated Estriol Reagent Wedge

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Symbol Definition

En: Upper limit of temperature ($\leq -20^{\circ}\text{C}$)
De: Obere Temperaturgrenze ($\leq -20^{\circ}\text{C}$)
Es: Límite superior de temperatura ($\leq -20^{\circ}\text{C}$)
Fr: Limite supérieure de température ($\leq -20^{\circ}\text{C}$)
It: Limite superiore di temperatura ($\leq -20^{\circ}\text{C}$)
Pt: Limite máximo de temperatura ($\leq -20^{\circ}\text{C}$)

2008-01

En: Lower limit of temperature ($\geq 2^{\circ}\text{C}$)
De: Mindesttemperatur ($\geq 2^{\circ}\text{C}$)
Es: Límite inferior de temperatura ($\geq 2^{\circ}\text{C}$)
Fr: Limite inférieure de température ($\geq 2^{\circ}\text{C}$)
It: Limite inferiore di temperatura ($\geq 2^{\circ}\text{C}$)
Pt: Limite mínimo de temperatura ($\geq 2^{\circ}\text{C}$)

En: Do not reuse
De: Nicht zur Wiederverwendung
Es: No reutilizar
Fr: Ne pas réutiliser
It: Non riutilizzare
Pt: Não reutilizar

En: Do not reuse
De: Nicht zur Wiederverwendung
Es: No reutilizar
Fr: Ne pas réutiliser
It: Non riutilizzare
Pt: Não reutilizar

En: Keep away from sunlight
De: Vor Sonneneinstrahlung schützen
Es: Proteger de la luz solar
Fr: Maintenir hors de portée de la lumière du soleil
It: Non esporre alla luce del sole
Pt: Manter afastado da luz solar

En: Batch code
De: Chargenbezeichnung
Es: Número de lote
Fr: Numéro de code de lot
It: Codice lotto
Pt: Código de lote

Symbol Definition

En: Contains sufficient for (n) tests
De: Es reicht für (n) Tests
Es: Contiene suficiente para (n) pruebas
Fr: Contient du matériel suffisant pour (n) tests
It: Contiene materiale sufficiente per (n) test
Pt: Contém o suficiente para (n) testes

En: Date format (year-month)
De: Datumformat (Jahr-Monat)
Es: Formato de fecha (año-mes)
Fr: Format de la date (année-mois)
It: Formato data (anno-mese)
Pt: Formato de data (ano-mês)

En: Use by
De: Verwendbar bis
Es: Fecha de caducidad
Fr: A utiliser avant
Pt: Usare entro
Pt: Usar até

En: Health Hazard
De: Gesundheitsgefährdung
Es: Peligro para la salud
Fr: Dangereux pour la santé
It: Pericolo per la salute
Pt: Perigo para a saúde

En: Exclamation Mark
De: Ausrufezeichen
Es: Signo de exclamación
Fr: Point d'exclamation
It: Punto esclamativo
Pt: Ponto de exclamação

En: Corrosion
De: Korrosion
Es: Corrosion
Fr: Corrosion
It: Corrosione
Pt: Corrosão

En: Skull and Crossbones
De: Totenkopf mit gekreuzten Knochen
Es: Calavera y tibias cruzadas
Fr: Tête de mort sur tibias croisées
It: Teschio e tibie incrociate
Pt: Caveira sobre tibias cruzadas

En: Environment
De: Umwelt
Es: Medio ambiente
Fr: Environnement
It: Ambiente
Pt: Ambiente

Limitations

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data representative of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:
ng/mL \times 3.467 \rightarrow nmol/L

Reportable Range: 0.07 – 12 ng/mL (0.24 – 42 nmol/L)

The assay calibration was determined based on the measurement of unconjugated estril in pregnancy serum (SST) samples by a gas chromatography-mass spectroscopy method (GC-MS) (see Method Comparison 1).

Analytical Sensitivity: Limit of Blank (highest value expected for a sample with no analyte; determined in accordance with CLSI EP17-A²⁶);
0.01 ng/mL (0.04 nmol/L)

Limit of Detection (lowest detectable concentration; determined in accordance with CLSI EP17-A²⁶);
0.04 ng/mL (0.14 nmol/L)

Functional Sensitivity: (concentration with 20% coefficient of variation (CV) determined in accordance with CLSI EP5-A²⁷);
0.15 ng/mL (0.52 nmol/L)

Gestational Week	Unconjugated Estril, ng/mL		
	n	Median	Central 95% Range
27	15	4.1	2.3 – 6.4
28	15	4.2	2.3 – 7.0
29	15	4.5	2.3 – 7.7
30	15	4.9	2.4 – 8.6
31	15	5.5	2.6 – 9.9
32	15	6.2	2.8 – 11.4
33	15	7.2	3.0 – >12
34	15	8.4	3.3 – >12
35	15	10.2	3.9 – >12
36	15	>12	4.7 – >12
37	15	>12	5.6 – >12
38	15	>12	6.6 – >12
39	15	>12	7.3 – >12
40	4	>12	7.6 – >12

Gestational Week	Unconjugated Estril, nmol/L		
	n	Median	Central 95% Range
27	15	14	8 – 22
28	15	15	8 – 24
29	15	16	8 – 27
30	15	17	8 – 30
31	15	19	9 – 34
32	15	21	10 – 40
33	15	25	10 – >42
34	15	29	11 – >42
35	15	35	14 – >42
36	15	>42	16 – >42
37	15	>42	19 – >42
38	15	>42	23 – >42
39	15	>42	25 – >42
40	4	>42	26 – >42

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.

© 2015 Siemens Healthcare Diagnostics. All rights reserved.

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PILL2KUE3 - 5

Changes in This Edition:

cd#EU22513: 1) Removed references to discontinued products. 2) Under Warnings and Precautions, added biohazard symbol and potential biohazard statement, caution statement, required GHS pictograms/signal words/hazard and precautionary codes and statements, offending chemicals and sources. 3) Under References, added citations from the potential biohazard statement. 4) Changed "Origin: UK" to "Made in: UK". 5) Updated the manufacturer address to include Glyn Rhonwy. 6) In Understanding the Symbols, revised the symbol for catalog number, added a GHS pictogram for health hazard, and updated existing symbols with GHS pictograms and terminology for exclamation mark, corrosion, skull and crossbones, and environment.

perinatal medicine. In: Klopfer A, editor. The human placenta. New York: Academic Press, 1980; 181-89. 19) Tulchinsky D. Use of biochemical indices in the management of high risk obstetric patients. Clin Perinatol 1980 Sep; 7(2):413-21. 20) Tulchinsky D. The role of maternal and fetal hormones in pregnancy. Contemp Obstet Gynecol 1981 Jun; 17(6):113-22. 21) Vorherr H. Placental insufficiency in relation to post-term pregnancy and fetal postmaturity. Amer J Obstet Gynecol 1975; 123:67-103. 22) Whipple MJ, et al. Estriol in pregnancy. VI: Experience with unconjugated plasma estriol assays and antepartum fetal heart rate testing in diabetic pregnancies. Amer J Obstet Gynecol 1979; 135:764-72. 23) DiCandolo JM. Urinary estriols and total serum estriol: a correlation study. Lab Med 1978; 9:47-8.

24) Goebelsmann U, et al. Estriol in pregnancy. Amer J Obstet Gynecol 1973; 115:795-802.

25) Townsley JD, et al. Circadian rhythms of serum and urinary estrogens in pregnancy. J Clin Endocrinol Metabol 1973; 36:289-95.

26) CLSI. Protocols for the Determination of Limits of Detection and Limits of Quantitation; Approved Guideline. CLSI document EP17-A Vol. 24 (No. 34). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2004. 27) Evaluation of Precision

Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition CLSI document EP5-A2. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2004. 28) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988; 37:377-82, 387-8.

29) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.

30) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor. www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Within-Run ¹			Total ²		
	Mean ³	SD ⁴	CV ⁵	SD	CV	CV
1	0.14	0.0209	14.9%	0.0275	19.6%	19.6%
2	0.32	0.0327	10.2%	0.0427	13.3%	13.3%
3	0.64	0.0510	8.0%	0.0680	10.6%	10.6%
4	1.17	0.0840	7.2%	0.1150	9.8%	9.8%
5	2.96	0.1690	5.7%	0.2580	8.7%	8.7%
6	5.97	0.3460	5.8%	0.4720	7.9%	7.9%
7	10.42	0.5420	5.2%	0.7610	7.3%	7.3%

Linearity (ng/mL)

	Dilution ¹	Observed ²		Expected ³		%O/E ⁴
		Mean	SD	Mean	SD	
1	8 in 8 ⁵	1.29	—	—	—	—
	7 in 8	1.21	1.13	1.13	107%	107%
	6 in 8	1.07	0.97	0.97	110%	110%
	5 in 8	0.90	0.81	0.81	111%	111%
	4 in 8	0.73	0.65	0.65	112%	112%
	3 in 8	0.55	0.48	0.48	115%	115%
	2 in 8	0.35	0.32	0.32	109%	109%
	1 in 8	0.17	0.16	0.16	106%	106%
2	8 in 8	2.43	—	—	—	—
	7 in 8	2.04	2.13	2.13	96%	96%
	6 in 8	1.80	1.82	1.82	99%	99%
	5 in 8	1.49	1.52	1.52	98%	98%
	4 in 8	1.26	1.22	1.22	103%	103%
	3 in 8	0.98	0.91	0.91	108%	108%
	2 in 8	0.63	0.61	0.61	103%	103%
	1 in 8	0.29	0.30	0.30	97%	97%
3	8 in 8	5.60	—	—	—	—
	7 in 8	4.44	4.90	4.90	91%	91%
	6 in 8	3.81	4.20	4.20	91%	91%
	5 in 8	3.35	3.50	3.50	96%	96%
	4 in 8	2.51	2.80	2.80	90%	90%
	3 in 8	2.05	2.10	2.10	98%	98%
	2 in 8	1.35	1.40	1.40	96%	96%
	1 in 8	0.69	0.70	0.70	99%	99%

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El Sistema de Calidad de Siemens Healthcare Diagnostics Products Ltd. está certificado por la ISO 13485:2003.